



913-906-9260

## DAILY FLUSHING OF PICC

### Gather Supplies

1 syringe Saline 4-5 cc

1 syringe Heparin 2-3 cc

alcohol swabs

1. Wash your hands.
2. Clean PICC line with x2 alcohol swabs-use friction and allow to dry
3. Attach syringe of Saline to PICC and flush using a gentle pulsing motion
4. Attach syringe of Heparin to end of PICC and flush using a gentle pulsing motion
5. Remove syringe and re-clamp PICC tubing (if clamp is present).

IF more than one PICC lumens: Flush **each** lumen of your PICC daily to ensure that it is still patent.