



913-906-9260

MEDICATION BY MEDISIS PUMP

1. Gather Supplies:

- Medisis pump
- 2 Syringes _____ Saline
- 1 Syringe _____ Heparin
- Tubing
- Alcohol wipes
- Medicine syringe

2. Preparing syringe for the pump:

Unscrew the pump all the way

Remove tubing from package (move white clamp towards center of tubing)

CLOSE white clamp

Connect white (female) end of tubing to colored end of syringe (must remove colored ends first)

Place syringe into Medisis pump (numbers and lines facing up)

Tighten the blue barrel **ALL THE WAY DOWN:** (you will feel resistance)

3. Priming the tubing:

Remove the white cap from the end of the tubing

Release the white clamp and watch for the fluid at the end of the tubing (this takes about 3-4 seconds)

Once you see the medication flow through to the end of the tubing, CLOSE white clamp

Place the protective white cap back on the end of the tubing

4. Starting the medication infusion:

Clean the end of the PICC line with 2 alcohol pads

Flush the PICC with _____cc saline

Remove the white cap from the medication tubing and connect to PICC line

Open the white clamp on the tubing –medicine will now infuse

5. Disconnecting after infusion is complete:

CLOSE the white clamp

Remove the tubing from the PICC line

Place small white cap on end of tubing if patient has more than one infusion per day of this specific medicine. (tubing can be reused for 24 hours if kept clean)

Flush PICC line with _____cc Saline

Flush PICC line with _____cc Heparin

—If using the tubing more than once in 24 hours place a new medicine syringe to the same tubing after the infusion is completed. You must remove the old syringe from the medisis first. Discard the old (now empty) syringe. Now place the tubing and new medicine syringe in the refrigerator for the next dose.

—Remove the medicine from the refrigerator 30 minutes before the dose is to be given

—Change the tubing every 24 hours